




INTRODUCING PERISOLV[®]: BIOFILM-ERASER



Elimination of the biofilm is key for the treatment success of periodontitis, peri-implant mucositis, and peri-implantitis.


BTC Speciality Health.

A NEW MINIMALLY INVASIVE APPLICATION

For periodontal and peri-implant inflammation

Periodontitis, peri-implant mucositis and peri-implantitis are bacterial inflammations with similar symptoms. The underlying cause of all three indications, which progress in a similar way, is bacterial plaque forming a biofilm, rich in pathogenic bacteria. The softening of the biofilm and effective elimination of the bacteria is thus the key prerequisite for effective treatment of these conditions.

Originated from Sweden, PERISOLV® is the first chemo mechanical treatment to be introduced to Australia. It is a new cleaning gel used in addition to mechanical debridement.



INDICATIONS:

PERISOLV® is indicated for use in:

- ✓ Peri-implantitis
- ✓ Mucositis
- ✓ Periodontitis

PERISOLV®-EFFECTS


- Elimination of the biofilm*
- Reduction of pocket depth even in persistent pockets^{1,2}

*Enhancing bacterial removal by mechanical debridement


Pocket depth reduction, even for persistent pockets^{1,2}

Positive in vitro data were also confirmed in a clinical case series study at the University of Ferrara.¹ The periodontal conditions of patients with infected deep residual pockets, a positive BOP and a pocket depth of ≥ 5 mm before treatment could be significantly improved after treatment with PERISOLV®. The average pocket depth before treatment was 5.7 ± 1.0 mm, decreasing to 3.4 ± 0.5 mm after treatment with ultrasound and PERISOLV®. After treatment, all pockets had a depth of ≤ 4 mm and were BOP-negative.

Average PPD



Average pocket depth**



**76% of 173 residual pockets

The effect of PERISOLV® on patients with residual pockets was documented via a case series in a private practice in Switzerland.² PERISOLV® was used in the treatment of 18 patients (173 residual pockets ≥ 5 mm). Follow-up ranged from 9 to over 30 weeks. A reduction in depth was determined in 76% of pockets. It was observed that in 66% of 5 mm pockets (depth from 5 mm to 3.87 mm), in 84% of 6 mm pockets (depth from 6 mm to 4.61 mm) and 92% of 7 mm pockets (depth from 7 mm to 5.62 mm), the probing depth was reduced. All other pockets remained stable. Using PERISOLV® also had a positive effect on BOP. Almost 70% of BOP-positive cases became BOP-negative.

MODE OF ACTION

Mechanical debridement alone or with PERISOLV®

Clear syringe with gel

- Water
- Carboxy methylcellulose (CMC)
- Sodium chloride
- Amino acids
- Titanium dioxide
- Alkaline pH

White syringe with liquid

- Water
- Sodium hypochlorite

PERISOLV®

TREATMENT OF PERI-IMPLANT MUCOSITIS

Case by Prof. Vincenzo Iorio-Siciliano, University of Catanzaro, Italy



Implant with probing depth (PD) $\leq$ 5mm and BOP⁺



Application of PERISOLV® before non-surgical treatment



Biofilm removal using a sonic scaler with PEEK tip



6 MONTHS AFTER TREATMENT
Probing depth (PD) after 6 months observation time

TREATMENT OF A FURCATION DEFECT

Case by Prof. Vincenzo-Iorio-Siciliano, University of Catanzaro, Italy




1. A PD of 5 mm was noted
2. A class-II furcation defect was recorded



3. First application of PERISOLV®
4. Scaling was performed using an ultrasonic device




5. Second application of PERISOLV®
6. Root planing was performed




- 6 MONTHS AFTER TREATMENT
7. A PD of 3 mm was reported
8. A class-I furcation defect with BOP was recorded


SIMPLE APPLICATION




1. CONNECT




2. MIX




3. CANNULA



4. APPLY



5. WAIT



6. START REGULAR SCALING





REFERENCES

1. Guarnelli ME et al. 'Professional local administration of chloramine-based treatment in conjunction with ultrasonic mechanical instrumentation: clinical outcomes in patients with deep periodontal pockets persisting following active non-surgical therapy'. Minerva Stomatologia, April 2015; Vol. 64 suppl. 1 al No. 2: 158-159.
2. Data on file 2016

PERISOLV® is a registered trademark and is manufactured by RLS Global AB, Neongatan 5, SE-431 53 Mölndal.

PERISOLV® is distributed by BTC Speciality Health Pty Ltd (a wholly owned company of BTC Health) under licence from RLS Global AB.

BTC Speciality Health.

Suite 201, 697 Burke Road, Camberwell VIC 3124, Australia.
Website: btchealth.com.au, Queries: queries@btchealth.com.au,
Phone No: 1800 100 282 PC: 0082018